

SAFETY DATA SHEET (EC 1907/2006)			
Ducera Liquid OCL universal			
Version:	2.3 / GB	Material no.	
Revision date:	24.11.2017	Specification	158262
Issue date:	07.11.2006	VA-Nr	01906947
replaces version:	2.2		
Page:	1 / 8		



SECTION 1: Identification of the substance/mixture and of the company/undertaking

1.1. Product identifier

Trade name	Ducera Liquid OCL universal
REACH Registration No.:	if available listed in Chapter. 3

1.2. Relevant identified uses of the substance or mixture and uses advised against

Relevant applications identified	For dental use only.
----------------------------------	----------------------

1.3. Details of the supplier of the safety data sheet

Company	DeguDent GmbH Postfach 1364 D-63403 Hanau
Telephone	+49 (0)6181/59-5576
Telefax	+49 (0)6181/59-5879
Email address	SDB.Degudent-DE@dentsplysirona.com

1.4. Emergency telephone number

Emergency information	+49 (0)6181/59-50 (This telephone number is available during office hours only.)
-----------------------	--

SECTION 2: Hazards identification

2.1. Classification of the substance or mixture

Classification according to Regulation (EC) No. 1272/2008 [CLP]

Not a hazardous mixture according to Regulation (EC) No 1272/2008.

2.2. Label elements

Labelling as per (EU) 1272/2008

Statutory basis	Labelling not required according to EU-CLP Ordinance (1272/2008).
-----------------	---

2.3. Other hazards

A PBT/vPvB evaluation is not available, since a chemical safety evaluation is not required / has not been carried out.

SECTION 3: Composition/information on ingredients

Chemical nature

The preparation contains:, aromatic oil

3.1. Substances

-

3.2. Mixtures

Information on ingredients / Hazardous components as per EU-CLP Regulation (EC) No. 1272/2008

• sodium hydroxide		< 0.15%	
CAS-No.	1310-73-2	EC-No.	215-185-5
Skin Corr. 1A			H314

SAFETY DATA SHEET (EC 1907/2006)

Ducera Liquid OCL universal

Version: **2.3 / GB**
Revision date: **24.11.2017**
Issue date: **07.11.2006**
replaces version: **2.2**
Page: **2 / 8**

Material no.
Specification **158262**
VA-Nr **01906947**



Texts of H phrases, see in Chapter 16

SECTION 4: First aid measures

4.1. Description of first aid measures

Remove contaminated or saturated clothing.

Inhalation

Move victims into fresh air.

Skin contact

Wash off with plenty of water.

Eye contact

With eye held open, thoroughly rinse immediately with plenty of water for at least 10 minutes.

Consult an ophthalmologist.

Ingestion

Have the mouth rinsed with water.

Obtain medical attention.

4.2. Most important symptoms and effects, both acute and delayed

Symptoms

No information available.

4.3. Indication of any immediate medical attention and special treatment needed

If skin sensitisation has developed and a causal relationship has been confirmed, further exposure should not be allowed.

SECTION 5: Firefighting measures

5.1. Extinguishing media

Suitable extinguishing media: Water spray
Foam
quenching powder
Carbon dioxide (CO₂)

Unsuitable extinguishing media: High volume water jet

5.2. Special hazards arising from the substance or mixture

The product itself does not burn.

5.3. Advice for firefighters

Employ protective equipment commonly used in the event of fire.

SECTION 6: Accidental release measures

6.1. Personal precautions, protective equipment and emergency procedures

Wear suitable gloves.

6.2. Environmental precautions

Do not allow the product into the following compartments: stretches of water, groundwater, sewage water

6.3. Methods and material for containment and cleaning up

Pick up mechanically with a suitable material and collect in a suitable container.

Disposal according to local authority regulations.

Suitable binder: universal absorbent

6.4. Reference to other sections

Wear personal protective equipment; see section 8.

Disposal considerations; see section 13.

SECTION 7: Handling and storage

7.1. Precautions for safe handling

SAFETY DATA SHEET (EC 1907/2006)

Ducera Liquid OCL universal

Version: **2.3 / GB**
Revision date: **24.11.2017**
Issue date: **07.11.2006**
replaces version: **2.2**
Page: **3 / 8**

Material no.
Specification **158262**
VA-Nr **01906947**



Local ventilation.
Wear suitable protective clothing and gloves.

7.2. Conditions for safe storage, including any incompatibilities

Advice on protection against fire and explosion

No special measures are required.

Storage

Keep tightly sealed in original packing.

Unsuitable materials Aluminium, tin, Zinc

German storage class

12 - Non Combustible Liquids

7.3. Specific end use(s)

We are unaware of any specific end uses which go beyond the data reported in Section 1.

SECTION 8: Exposure controls/personal protection

8.1. Control parameters

• sodium hydroxide

CAS-No.	1310-73-2	EC-No.	215-185-5
Control parameters	2 mg/m3		Short Term Exposure Limit (STEL):(EH40 WEL)

8.2. Exposure controls

Engineering measures

Local ventilation.

Personal protective equipment

Hand protection

Preventive skin protection, Use barrier cream regularly.

The suitability for a specific workplace should be discussed with the producers of the protective gloves.,
The exact break through time can be obtained from the protective glove producer and this has to be observed.

Glove material Fluorinated rubber (FKM)

Material thickness 0.4 mm

Break through time 480 min

Method Source: GESTIS substance database (hazardous substance information system of commercial professional associations)

Glove material butyl-rubber

Material thickness 0.5 mm

Break through time 480 min

Method Source: GESTIS substance database (hazardous substance information system of commercial professional associations)

Wear protective gloves made of the following materials:., alkali-resistant

Glove material Polychloroprene (PCP)

Material thickness 0.5 mm

Break through time 480 min

Method Source: GESTIS substance database (hazardous substance information system of commercial professional associations)

Glove material PVC

Material thickness 0.5 mm

Break through time 480 min

Method Source: GESTIS substance database (hazardous substance information system of commercial professional associations)

Eye/face protection

Safety glasses with side-shields

SAFETY DATA SHEET (EC 1907/2006)**Ducera Liquid OCL universal**

Version: **2.3 / GB**
Revision date: **24.11.2017**
Issue date: **07.11.2006**
replaces version: **2.2**
Page: **4 / 8**

Material no.
Specification **158262**
VA-Nr **01906947**

**Hygiene measures**

Handle in accordance with good industrial hygiene and safety practice., Do not eat, drink, smoke, or sniff while at work. Wash your hands and/or face before breaks and before termination of work.

Protective measures

If there is the possibility of skin/eye contact, the indicated hand/eye/body protection should be used.

SECTION 9: Physical and chemical properties**9.1. Information on basic physical and chemical properties****Appearance**

Form liquid
Colour colourless

Odour fruity

Odour threshold: no data available

pH ca. 13 Medium: Water

Melting point/range ca. 0 °C

Boiling point/range ca. 100 °C

Flash point not applicable

Evaporation rate no data available

Flammability (solid, gas) not applicable

Lower explosion limit not applicable

Upper explosion limit not applicable

Density ca. 1.1 g/cm³

Water solubility soluble

Partition coefficient: n-octanol/water no data available

Autoinflammability Not capable of spontaneous combustion or heating.

Thermal decomposition no data available

Viscosity, dynamic no data available

Explosiveness not explosive

Oxidizing properties not oxidizing

9.2. Other information

Metal corrosion Not expected during handling from practical experience.

Other information No further physicochemical data were determined.

SECTION 10: Stability and reactivity**10.1. Reactivity**

no data available

SAFETY DATA SHEET (EC 1907/2006)				
Ducera Liquid OCL universal				
Version:	2.3 / GB	Material no.		
Revision date:	24.11.2017	Specification	158262	
Issue date:	07.11.2006	VA-Nr	01906947	
replaces version:	2.2			
Page:	5 / 8			

- 10.2.
Chemical stability

The product is chemically stable.
- 10.3.
Possibility of hazardous reactions

Possibility of hazardous reactions

No dangerous reactions known.
- 10.4.
Conditions to avoid

None known.
- 10.5.
Incompatible materials

Acids, Ammonium salts
- 10.6.
Hazardous decomposition products

None known

SECTION 11: Toxicological information

11.1. Information on toxicological effects

Animal experiment studies focusing on this formulation are not available.

Acute oral toxicity	no data available
Acute inhalation toxicity	no data available
Acute dermal toxicity	no data available
Skin irritation	no data available
Eye irritation	no data available
Sensitization	no data available
Repeated dose toxicity	no data available
Mutagenicity assessment	no data available
Carcinogenicity	No data available
Toxicity to reproduction	No data available

SECTION 12: Ecological information

12.1. Toxicity

No ecotoxicological data is available for this product.

12.2. Persistence and degradability

Biodegradability	no data available
------------------	-------------------

12.3. Bioaccumulative potential

Bioaccumulation	no data available
-----------------	-------------------

12.4. Mobility in soil

Mobility	No data available
----------	-------------------

SAFETY DATA SHEET (EC 1907/2006)**Ducera Liquid OCL universal**

Version: **2.3 / GB**
Revision date: **24.11.2017**
Issue date: **07.11.2006**
replaces version: **2.2**
Page: **6 / 8**

Material no.
Specification **158262**
VA-Nr **01906947**

**12.5. Results of PBT and vPvB assessment**

A PBT/vPvB evaluation is not available, since a chemical safety evaluation is not required / has not been carried out.

12.6. Other adverse effects

Further Information

Dusts and water-soluble forms of the alloy:, Introduction into soil, natural water bodies or sewerage must be prevented.

SECTION 13: Disposal considerations**13.1. Waste treatment methods****Product**

Disposal according to local authority regulations.

Uncleaned packaging

Disposal according to local authority regulations.

SECTION 14: Transport information

Not dangerous according to transport regulations.

- | | |
|-------------------------------------|----|
| 14.1. UN number: | -- |
| 14.2. UN proper shipping name: | -- |
| 14.3. Transport hazard class(es): | -- |
| 14.4. Packing group: | -- |
| 14.5. Environmental hazards: | -- |
| 14.6. Special precautions for user: | No |

SECTION 15: Regulatory information**15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture****National legislation****15.2. Chemical safety assessment**

Chemical safety assessment No Chemical Safety Report as per Articles 2(8), 2(9) or 14 of the REACH Regulation is required for this product.

SECTION 16: Other information

Classification and applied procedure to derive the classification of mixtures according to EU Regulation (EC) No. 1272/2008 (CLP)

Relevant H phrases from chapter 3

H314 : Causes severe skin burns and eye damage.

Further information

Changes since the last version are highlighted in the margin. This version replaces all previous versions.

SAFETY DATA SHEET (EC 1907/2006)

Ducera Liquid OCL universal

Version: **2.3 / GB**
Revision date: **24.11.2017**
Issue date: **07.11.2006**
replaces version: **2.2**
Page: **7 / 8**

Material no.
Specification **158262**
VA-Nr **01906947**



This information and all further technical advice is based on our present knowledge and experience. However, it implies no liability or other legal responsibility on our part, including with regard to existing third party intellectual property rights, especially patent rights. In particular, no warranty, whether express or implied, or guarantee of product properties in the legal sense is intended or implied. We reserve the right to make any changes according to technological progress or further developments. The customer is not released from the obligation to conduct careful inspection and testing of incoming goods. Performance of the product described herein should be verified by testing, which should be carried out only by qualified experts in the sole responsibility of a customer. Reference to trade names used by other companies is neither a recommendation, nor does it imply that similar products could not be used.

Legend

ADR	European Agreement concerning the International Carriage of Dangerous Goods by Road
ADN	European Agreement concerning the International Carriage of Dangerous Goods by Inland Waterways
ASTM	American Society for Testing and Materials
ATP	Adaptation to Technical Progress
BCF	Bioconcentration factor
BetrSichV	German Ordinance on Industrial Safety and Health
c.c.	closed cup
CAS	Chemical Abstract Services
CESIO	European Committee of Organic Surfactants and their Intermediates
ChemG	German Chemicals Act
CMR	carcinogenic-mutagenic-toxic for reproduction
DIN	German Institute for Standardization
DMEL	Derived minimum effect level
DNEL	Derived no effect level
EINECS	European Inventory of Existing Commercial Chemical Substances
EC50	half maximal effective concentration
GefStoffV	German Ordinance on Hazardous Substances
GGVSEB	German ordinance for road, rail and inland waterway transportation of dangerous goods
GGVSee	German ordinance for sea transportation of dangerous goods
GLP	Good Laboratory Practice
GMO	Genetic Modified Organism
IATA	International Air Transport Association
ICAO	International Civil Aviation Organization
IMDG	International Maritime Dangerous Goods
ISO	International Organization For Standardization
LOAEL	Lowest observed adverse effect level
LOEL	Lowest observed effect level
NOAEL	No observed adverse effect level
NOEC	no observed effect concentration
NOEL	no observed effect level
o. c.	open cup
OECD	Organisation for Economic Cooperation and Development
OEL	Occupational Exposure Limit
PBT	Persistent, bioaccumulative, toxic
PEC	Predicted effect concentration
PNEC	Predicted no effect concentration
REACH	REACH registration
RID	Convention concerning International Carriage by Rail
STOT	Specific Target Organ Toxicity
SVHC	Substances of Very High Concern
TA	Technical Instructions
TPR	Third Party Representative (Art. 4)
TRGS	Technical Rules for Hazardous Substances
VCI	German chemical industry association
vPvB	very persistent, very bioaccumulative

SAFETY DATA SHEET (EC 1907/2006)				
Ducera Liquid OCL universal				
Version:	2.3 / GB	Material no.		
Revision date:	24.11.2017	Specification	158262	
Issue date:	07.11.2006	VA-Nr	01906947	
replaces version:	2.2			
Page:	8 / 8			

VOC	volatile organic compounds
VwVwS	German Administrative Regulation on the Classification of Substances Hazardous to Waters into Water Hazard Classes
WGK	Water Hazard Class
WHO	World Health Organization