

Trade name: **TR-Neutralizer** Issue/Revision: 07.01.2021 Print date: 15.01.2021

Version: 2.0 Replaces version: 1.0

### SECTION 1: Identification of the substance/mixture and of the company/undertaking 1.1. Product identifier Trade name: **TR-Neutralizer** 1.2. Relevant identified uses of the substance or mixture and uses advised against Relevant identified uses: Neutralizing and cleaning agent Intended purpose: Liquid concentrate for the neutralization after the alkaline cleaning in washer-disinfectors. Uses advised against: None at intended use. Note: The product is intended for professional users. 1.3. Details of the supplier of the safety data sheet Manufacturer/Supplier: ALPRO MEDICAL GMBH Mooswiesenstraße 9 D-78112 St. Georgen (Germany) Telephone: +49 7725 9392-0 Telefax: +49 7725 9392-91 E-mail: alpro@alpro-medical.de Internet: www.alpro-medical.com E-mail address for the competent person responsible for the safety data sheet: doku@alpro-medical.de 1.4. Emergency telephone number In-house emergency telephone number: +49 7725 9392-0 Monday – Friday from 08:00 am to 04:30 pm (UTC+1); for chemical information and legal information on hazardous substances only +49 761 19240 Poison centre: Poisoning information centre, Freiburg, Germany (24 h / 7 d), English is spoken

# SECTION 2: Hazards identification

## 2.1. Classification of the substance or mixture

## Classification in accordance with Regulation (EC) No 1272/2008 [CLP]

Classification	Classification procedure
Skin Corr. 1B; H314	Calculation method

Full text of hazard classes as well as H-phrases: see under SECTION 16.1.

### 2.2. Label elements

## Label elements in accordance with Regulation (EC) No 1272/2008 [CLP]

Hazard pictograms:



Signal word:



Trade name: **TR-Neutralizer** 

Issue/Revision: 07.01.2021 Print date: 15.01.2021

Hazard compo for labelling:	nents	50 - 70 % phosphoric acid (7664-38-2)
H-phrases:	H314	Causes severe skin burns and eye damage.
P-phrases:	P280	Wear protective gloves/protective clothing/eye protection/face protection.
	P301+P330+P331	IF SWALLOWED: rinse mouth. Do NOT induce vomiting.
	P303+P361+P353	IF ON SKIN (or hair): Remove/Take off immediately all
		contaminated clothing. Rinse skin with water/shower.
	P305+P351+P338	IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
	P310	Immediately call a POISON CENTER/ doctor.

## 2.3. Other hazards

The substances in the mixture do not meet the PBT/vPvB criteria according to Regulation (EC) No 1907/2006, Annex XIII.

No further hazards known.

# **SECTION 3: Composition/information on ingredients**

## 3.2. Mixtures

Chemical characterisation:

Mixture of substances listed below with non-hazardous additions in aqueous solution.

### Hazardous ingredients

Chemical name	-	Classification in accordance with Regulation (EC) No 1272/2008	Weight %
Phosphoric acid	CAS No: 7664-38-2 EC No: 231-633-2 Index No: 015-011-00-6 REACH Registration No: 01-2119485924-24-XXXX	Skin Corr. 1B; H314 Skin Irrit. 2; H315 Eye Irrit. 2; H319 Specific concentration limits: Skin Corr. 1B; H314: $C \ge 25 \%$ Skin Irrit. 2; H315: 10 $\% \le C < 25 \%$ Eye Irrit. 2; H319: 10 $\% \le C < 25 \%$	≥ 50 - < 70

Full text of hazard classes and H-phrases: see SECTION 16.1.

# **SECTION 4: First aid measures**

## 4.1. Description of first aid measures

General information:	First aider: Pay attention to self-protection! Remove contaminated, saturated clothing immediately.
Following inhalation:	Move affected person into fresh air and keep still and warm. Seek medical advice.
Following skin contact:	Wash skin immediately with plenty of water and soap. In case of skin reactions, consult a physician.
Following eye contact:	Flush eyes immediately with flowing water for 10 to 15 minutes holding eyelids apart. Remove contact lenses, if present and easy to do. Consult an ophthalmologist.



Replaces version: 1.0

Version: 2.0

Trade name: **TR-Neutralizer** Issue/Revision: 07.01.2021

Print date: 15.01.2021

Following ingestion: Rinse mouth with water. Let drink plenty of water. Do not induce vomiting (risk of perforation). Consult a physician immediately.

### 4.2. Most important symptoms and effects, both acute and delayed

Causes severe skin burns and eye damage.

## 4.3. Indication of any immediate medical attention and special treatment needed

No information available.

### **SECTION 5: Firefighting measures**

### 5.1. Extinguishing media

Suitable extinguishing media: Water spray jet, alcohol resistant foam, extinguishing powder, carbon dioxide (CO<sub>2</sub>)

Unsuitable extinguishing media: Full water jet

### 5.2. Special hazards arising from the substance or mixture

Hazardous combustion products: Phosphorus oxides (P<sub>x</sub>O<sub>y</sub>)

## 5.3. Advice for firefighters

Special protective equipment:Wear self-contained breathing apparatus.Further information:Cool endangered containers with water spray jet. Collect<br/>contaminated fire extinguishing water separately. Do not discharge<br/>into drains or rivers.

# **SECTION 6: Accidental release measures**

### 6.1. Personal precautions, protective equipment and emergency procedures

### For non-emergency personnel

Use personal protective equipment. See SECTION 8.2.

Avoid skin and eye contact. Special danger of slipping by leaked/spilled product. Evacuate danger area. Observe emergency plans. Consult experts.

### For emergency responders

Use personal protective equipment. See SECTION 8.2.

### 6.2. Environmental precautions

Do not discharge into drains or rivers.

### 6.3. Methods and material for containment and cleaning up

### Containment

For large spills, dyke spilled material or otherwise contain material to ensure runoff does not reach a waterway. Cover or seal drains.

### Cleaning up

Wipe up small amounts with absorbent material (e.g. cloth, fleece). Absorb large amounts with liquidbinding material (sand, diatomaceous earth, universal binder, sawdust). Collect in suitable, closed containers for disposal. Clean contaminated surfaces thoroughly.

### Other information

Inappropriate containment and cleaning methods are not known.



Trade name: **TR-Neutralizer** 

Issue/Revision: 07.01.2021 Print date: 15.01.2021 Version: 2.0 Replaces version: 1.0

### 6.4. Reference to other sections

Information on safe handling see SECTION 7.1. Information on personal protective equipment see SECTION 8.2. Information on disposal see SECTION 13.

### **SECTION 7: Handling and storage**

### 7.1. Precautions for safe handling

### Precautions

Avoid contact with skin and eyes.

### Advice on general occupational hygiene

When using do not eat, drink or smoke. Wash hands before breaks and at end of work. Keep away from food and drink.

### 7.2. Conditions for safe storage, including any incompatibilities

Requirements for storage rooms and vessels:	Keep only in the original container. Keep container tightly closed and kept upright to prevent any leakage.
Advice on common storage:	Keep away from food, drink and animal feedingstuffs.
Further information on storage conditions:	Not necessary
Storage class ([DE] TRGS 510):	LGK 8B Non-combustible corrosive hazardous substances

### 7.3. Specific end use(s)

Apart from the uses mentioned in section 1.2 no other specific end uses are stipulated.

### Industry and sector specific guidance

- [DE] TRGS 525 Hazardous substances in medical care facilities (Section 7 Activities with disinfectants); Issue: September 2014; Source: GMBI 2014 page 1294-1307 of 13.10.2014 [No 63]; www.baua.de
- [DE] DGUV rules 107-002 (former BGR 206) Disinfection works in health service Issue: July 1999; Source: www.dguv.de/publikationen

### **SECTION 8: Exposure controls/personal protection**

### 8.1. Control parameters

### **Occupational exposure limits**

Limit values						
Country	Long term (8 hours)		Short term (15 minutes)		Legal basis	Remarks
	ррт	mg/m³	ррт	mg/m³		
Phosphoric acid (CAS No: 7664-38-2)						
EU	-	1	-	2	2000/39/EC	-
UK	-	1	-	2	EH40	

Used abbreviations, symbols, numerals and explanations in column "Remarks"

### **Biological limit values**

Does not contain substances above concentration limits fixing a biological limit value.



# Trade name: **TR-Neutralizer**

Issue/Revision: 07.01.2021 Print date: 15.01.2021

Version: 2.0 Replaces version: 1.0

### Information on monitoring procedures

BS EN 482:2012-04-30; Title: Workplace exposure - General requirements for the performance of procedures for the measurement of chemical agents; British version of EN 482:2012

BS EN 689:1996-04-15; Title: Workplace atmospheres - Guidance for the assessment of exposure by inhalation to chemical agents for comparison with limit values and measurement strategy; British version of EN 689:1995

BS EN 14042:2003-04-24; Title: Workplace atmospheres - Guide for the application and use of procedures for the assessment of exposure to chemical and biological agents; British version of EN 14042:2003

### 8.2. Exposure controls

### Appropriate engineering controls

### Technical and organisational protective measures

The eyewash station (or eyewash bottle) and emergency shower must be located near the workplace.

# Personal protective equipment

Eye/face protection:	Safety glasses with side protection according to EN 166
Skin protection:	
Hand protection:	Protective gloves according to EN 374 <u>Splash guard:</u> Disposable gloves made of nitrile rubber (thickness 0.11 mm) <u>Permanent contact (&gt; 480 min):</u> Protective gloves made of nitrile rubber (thickness 0.40 mm)
Other skin protection:	Long-sleeved protective clothing (lab coat)
Respiratory protection:	Not necessary when used as intended.
Thermal hazards:	No special protective measures necessary.

### **Environmental exposure controls**

Do not discharge into drains.

### **SECTION 9: Physical and chemical properties**

### 9.1. Information on basic physical and chemical properties

Appearance:	clear, colourless liquid	
Odour:	odourless	
Odour threshold:	no data available	
pH (undiluted):	< 1	(20 °C)
Melting point/freezing point:	no data available	
Initial boiling point and boiling range:	no data available	
Flash point:	not applicable	
Evaporation rate:	no data available	
Flammability (solid, gas):	not applicable	
Lower explosive limit:	not applicable	
Upper explosive limit:	not applicable	

# Trade name: **TR-Neutralizer**

Issue/Revision: 07.01.2021 Print date: 15.01.2021



Version: 2.0 Replaces version: 1.0

	Vapour pressure:	no data available	( °C)
	Vapour density:	no data available	
	Relative density:	1.35 - 1.45	(20 °C)
	Solubility in water:	completely soluble	
	Partition coefficient: n-octanol/water	not applicable	
	Auto-ignition temperature:	not applicable	
	Decomposition temperature:	no data available	
	Viscosity:	no data available	
	Explosive properties:	none	
	Oxidising properties:	none	
9.2.	Other information		
	Refractive index nD:	no data available	( °C)
	Electrical conductivity (undiluted):	3500 - 4000 μS/cm	(20 °C)

# **SECTION 10: Stability and reactivity**

### 10.1. Reactivity

No hazardous reactions when handled and stored as intended.

### 10.2. Chemical stability

The product is stable when handled and stored as intended.

## 10.3. Possibility of hazardous reactions

None known

### 10.4. Conditions to avoid

None known

# 10.5. Incompatible materials

None known

## 10.6. Hazardous decomposition products

Does not decompose when used as intended.

# **SECTION 11: Toxicological information**

# 11.1. Information on toxicological effects

# Acute toxicity

### Product

Acute toxicity - oral:	Acute Toxicity Estimate ATE <sub>mix</sub> > 2000 mg/kg => no classification
Acute toxicity - dermal:	Acute Toxicity Estimate ATE <sub>mix</sub> > 2000 mg/kg => no classification
Acute toxicity - inhalation:	Acute Toxicity Estimate ATE <sub>mix</sub> > 20 mg/l => no classification

# Trade name: **TR-Neutralizer**

Issue/Revision: 07.01.2021 Print date: 15.01.2021



Version: 2.0 Replaces version: 1.0

# Skin corrosion/irritation

Product

Causes severe skin burns. [calculation method]

### Serious eye damage/irritation

### Product

Causes serious eye damage. [calculation method]

### Respiratory or skin sensitisation

Product

No data available.

### Germ cell mutagenicity

Product

No data available.

### Carcinogenicity

Product

No data available.

### **Reproductive toxicity**

Product

No data available.

### **STOT-single exposure**

Product

No data available.

# STOT-repeated exposure

Product

No data available.

### Aspiration hazard

Product

No data available.

## **SECTION 12: Ecological information**

### 12.1. Toxicity

No data available.

### 12.2. Persistence and degradability

Biodegradability:

Phosphoric acid as inorganic substance is not subject to biodegradation. At normal dilution it does not endanger the function of a biological wastewater treatment plant and it will be eliminated for example in a phosphate precipitation stage.

### 12.3. Bioaccumulative potential

No data available.



# Trade name: **TR-Neutralizer**

Issue/Revision: 07.01.2021 Print date: 15.01.2021

Version: 2.0 Replaces version: 1.0

### 12.4. Mobility in soil

No data available.

## 12.5. Results of PBT and vPvB assessment

The substances in the mixture do not meet the PBT/vPvB criteria according to Regulation (EC) No 1907/2006, Annex XIII.

### 12.6. Other adverse effects

No data available.

### **SECTION 13: Disposal considerations**

#### 13.1. Waste treatment methods

### **Disposal of the product**

Product residues must be disposed of as hazardous waste in compliance with the Directive 2008/98/EC on waste as well as national and regional regulations. Do not dispose of via the waste water. Leave product in the original container as possible. Do not mix with other waste materials.

Waste codes / waste designations according to EWC

Product residues: 06 01 04\* phosphoric and phosphorous acid

#### Disposal of the packaging

Packaging contaminated with product is considered as hazardous waste and must be disposed of accordingly.

Waste codes / waste designations according to EWC

Contaminated packaging:	15 01 10* packaging containing residues of or contaminated by
	hazardous substances

#### Recommendation

Contaminated packaging must be emptied optimally and can be recycled after appropriate cleaning (rinse with water).

### **SECTION 14: Transport information**

### 14.0. Transport classification

Dangerous good in sense of the transport regulations in road traffic (ADR), railway traffic (RID), inland waterway traffic (ADN), maritime traffic (IMDG-Code) and air traffic (ICAO-TI/IATA-DGR).

#### 14.1.UN number

UN 1805

### 14.2. UN proper shipping name

### ADR/RID/ADN

PHOSPHORIC ACID, SOLUTION

### IMDG-Code/ICAO-TI/IATA-DGR

### PHOSPHORIC ACID, SOLUTION

# 14.3. Transport hazard class(es)

Class: 8 Subsidiary risk(s): -

Safety data sheet according to Regulation (EC) No 1907/2006	
Trade name: TR-Neutralizer	
Issue/Revision: 07.01.2021 Print date: 15.01.2021	Version: 2.0 Replaces version: 1.0
14.4. Packing group	
14.5. Environmental hazards	
ADR/RID/ADN	
Environmentally Hazardous:	No
IMDG-Code	
Marine Pollutant:	No
14.6. Special precautions for user	
Not necessary.	
14.7. Transport in bulk according to Annex II of MARPOL 7	3/78 and the IBC Code
Not applicable for product as supplied.	
14.8. Further information	
Transport category according to ADR section 1.1.3.6:	3
Maximum total quantity per transport unit according to ADR section 1.1.3.6:	1000 L
Limited quantity (Maximum quantity per inner packaging) according to ADR/RID/ADN/IMDG-Code:	5 L
Classification code according to ADR/RID/ADN:	C1
Hazard identification number according to ADR/RID:	80
Tunnel restriction code according to ADR/RID:	E
Segregation group according to IMDG-Code section 5.4.1.5.11.1:	IMDG-Code- Segregation group 1 – acids
EmS codes:	F-A, S-B

# **SECTION 15: Regulatory information**

### 15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture

# **EU-Regulations**

REGULATION (EC) No 1005/2009 on substances that deplete the ozone layer

not applicable

REGULATION (EC) No 850/2004 on persistent organic pollutants and amending Directive 79/117/EEC not applicable

REGULATION (EU) No 649/2012 concerning the export and import of hazardous chemicals

not applicable

DIRECTIVE 2012/18/EU (Seveso III Directive) on the control of major-accident hazards involving dangerous substances, amending and subsequently repealing Council Directive 96/82/EC

not applicable

DIRECTIVE 2010/75/EU on industrial emissions (integrated pollution prevention and control)

not applicable



# Trade name: **TR-Neutralizer**

Issue/Revision: 07.01.2021 Print date: 15.01.2021 Version: 2.0 Replaces version: 1.0

REACH - List of substances subject to authorisation (Annex XIV)

not applicable

REACH – Restrictions on the manufacture, placing on the market and use of certain dangerous substances, mixtures and articles (Annex XVII)

not applicable

COUNCIL DIRECTIVE 94/33/EC on the protection of young people at work

Observe employment restrictions for juveniles.

COUNCIL DIRECTIVE 92/85/EEC on the introduction of measures to encourage improvements in the safety and health at work of pregnant workers and workers who have recently given birth or are breastfeeding

not applicable

### 15.2. Chemical safety assessment

For this mixture no chemical safety assessment has been carried out.

## **SECTION 16: Other information**

### 16.1. Full text of hazard classes and H-phrases

### Hazard classes

Eye Irrit.	Eye irritation
Skin Corr.	Skin corrosion
Skin Irrit.	Skin irritation

### H-phrases (Hazard statements)

- H314Causes severe skin burns and eye damage.H315Causes skin irritation.
- H319 Causes serious eye irritation.

### 16.2. Abbreviations and acronyms

- ADN <u>Accord européen relatif au transport international des marchandises d</u>angereuses par voie de <u>n</u>avigation intérieure (European Agreement concerning the International Carriage of Dangerous Goods by Inland Waterways)
- ADR <u>Accord européen relatif au transport international des marchandises d</u>angereuses par <u>r</u>oute (European Agreement concerning the International Carriage of Dangerous Goods by Road)
- BGRBerufsgenossenschaftliche Regeln (English: Employers' liability insurance association rules)BSBritish Standards
- CAS <u>Chemical Abstracts Service</u>
- CLP Regulation on <u>Classification</u>, <u>Labelling and Packaging of Substances and Mixtures</u>
- [DE] National German regulations
- DGUV <u>D</u>eutsche <u>G</u>esetzliche <u>U</u>nfall<u>v</u>ersicherung (English: German statutory accident insurance) EC European Community
- EEC <u>European Economic Community</u>
- EmS <u>Emergency Schedules (Emergency response procedures for ships carrying dangerous goods)</u>
- EN European Standard
- EU <u>European Union</u>
- EWC <u>European Waste Catalogue</u>



Trade name: **TR-Neutralizer** 

Issue/Revision: 07.01.2021 Print date: 15.01.2021 Version: 2.0 Replaces version: 1.0

	GHS	<u>G</u> lobally <u>H</u> armonized <u>System of Classification</u> , Labelling and Packaging of Chemicals
	GMBI	<u>G</u> emeinsames <u>M</u> inisterial <u>bl</u> att (English: Joint Ministerial Gazette)
	IATA-DGR	International <u>Air Transport Association - Dangerous Goods Regulations</u>
	IBC-Code	International Code for the Construction and Equipment of Ships carrying Dangerous Chemicals in Bulk
	ICAO-TI	Technical Instructions For The Safe Transport of Dangerous Goods by Air
	IMDG-Code	International Maritime Code for Dangerous Goods
	LGK	Lagerklasse (English: Storage class)
	MARPOL	International Convention for the Prevention of Marine Pollution from Ships
	PBT	Persistent, <u>b</u> ioaccumulative and <u>t</u> oxic
	ppm	<u>P</u> arts <u>p</u> er <u>m</u> illion
	REACH	Registration, Evaluation, Authorisation and Restriction of Chemicals
	RID	<u>R</u> èglement concernant le transport <u>International ferroviaire de marchandises</u> <u>D</u> angereuses
		(Regulations Concerning the International Carriage of Dangerous Goods by Rail)
	TRGS	<u>Technische Regeln für Gefahrstoffe (English: Technical Rules for Hazardous Substances)</u>
	UN	<u>U</u> nited <u>N</u> ations
	UTC	Coordinated Universal Time (French: Temps Universel Coordonné)
	vPvB	<u>V</u> ery <u>p</u> ersistent and <u>v</u> ery <u>b</u> ioaccumulative
16.3. Key literature references and sources for data		
- Regulation (EC) No 1907/2006 (REACH), Annex II		

- European Chemicals Agency (ECHA) Guidance on the compilation of safety data sheets; Version 2.1 (February 2014); http://echa.europa.eu/documents/10162/13643/sds\_en.pdf
- GISBAU (Hazardous substances information system of the BG BAU) course "safety data sheet"; <a href="http://www.bgbau.de/gisbau/SDB/lehrgang/lehrgang.htm">http://www.bgbau.de/gisbau/SDB/lehrgang/lehrgang.htm</a>
- Regulation (EC) No 1272/2008 (CLP regulation)
- European Chemicals Agency (ECHA) Guidance on Labelling and Packaging in accordance with Regulation (EC) No 1272/2008 (04/2011);
  - http://echa.europa.eu/documents/10162/13562/clp\_labelling\_en.pdf
- European Chemicals Agency (ECHA), Registered substances; <u>http://echa.europa.eu/information-on-chemicals/registered-substances</u>
- European Chemicals Agency (ECHA), C&L Classification and Labelling Inventory; <u>http://echa.europa.eu/information-on-chemicals/cl-inventory-database</u>
- Institute for Occupational Safety and Health of the German Social Accident Insurance (IFA):
  GESTIS database on hazardous substances and GESTIS International limit values for chemical agents;
  <a href="http://www.dguv.de/dguv/ifa/index.jsp">http://www.dguv.de/dguv/ifa/index.jsp</a>
- German Environmental Agency (Umweltbundesamt), Section IV 2.4: Office of Documentation and Information on Substances Hazardous to Waters RIGOLETTO (catalogue of Substances Hazardous to Waters); <u>http://webrigoletto.uba.de/rigoletto</u>

# 16.4. Training advice

Provide adequate information, instructions and training for users.

# 16.5. Indication of changes

A dash in the left hand margin indicates an amendment from the previous version.

The information given in the safety data sheet only applies to the described product in connection with its intended use. This information is based on the latest state of our knowledge at the time of revision. In particular, it describes our product under the aspect of its hazards and safety measures to be taken. It does not constitute any guarantee of product properties and quality features.