

page 1/6 BeautiSealant PASTE

Printing date: June 10, 2020

SECTION 1. Identification of the substance or mixture and of the supplier

1.1 Product identifier

Trade Name:

BeautiSealant "PASTE"

1.2 Relevant identified uses of the substance or mixture and uses advised against

Relevant identified uses: Dental material Uses advised against: No further data

1.3 Details of the supplier of the safety date sheet

Company/Undertaking identification

Manufacturer's Name: SHOFU DENTAL GmbH Address: An der Pönt 70, 40885 Ratingen, Germany

Phone: +49 (0) 2102-8664-0 Fax: +49 (0) 2102-8664-64 E-Mail: info@shofu.de

Section in Charge: Quality Management & Regulatory Affairs

1.4 Emergency Telephone Number

+49-2102-8664-53 (SHOFU DENTAL GmbH) 24 hours / 7 days

SECTION 2. Hazards identification

2.1 Classification of the substance or mixture

CLASSIFICATION (EC 1272/2008)

Skin Irrit. 2 H315 Causes skin irritation

Skin Sens. 1 H317 May cause an allergic skin reaction

Eye Irrit. 2 H319 Causes serious eye irritation

2.2 Label elements

LABEL IN ACCORDANCE WITH (EC) NO.1272/2008



GHS07

HAZARD-DETERMINING COMPONENTS OF LABELLING

UDMA

SIGNAL WORD

Warning

HAZARD STATEMENTS

H315 Causes skin irritation.

H317 May cause an allergic skin reaction.

H319 Causes serious eye irritation.

PRECAUTIONARY STATEMENTS

P280 Wear protective gloves/protective clothing/eye protection/face

protection.

(Contd. on page2)

Version Number 5

Revision date: April 1, 2020



page 2/6 BeautiSealant PASTE

Printing date: June 10, 2020

(Contd. of page1)

P302+P352

IF ON SKIN: Wash with plenty of soap and water.

P305+P351+P338

IF IN EYES: Rinse cautiously with water for several minutes.

Remove contact lenses, if present and easy to do. Continue rinsing.

P501 Dispose of contents/container in accordance with

local/regional/national/international regulation.

2.3 Other hazards

Results of PBT and vPvB assessment

PBT: Not applicable. vPvB: Not applicable.

SECTION 3. Composition/information on ingredients

3.1 Chemical characterization: Mixtures

3.2 Description: Mixture of substances listed below with nonhazardous additions.

3.3 Dangerous components:

3.4 Additional information: For the wording of the listed risk phrases refer to section 2.

SECTION 4. First-aid measures

4.1 Description of first aid measures

Eye contact: Rinse cautiously with water for several minutes. Remove contact lenses, if

present. and easy to do. If eye irritation persists, get medical advice/attention.

Skin contact: Wash immediately with soap and plenty of water. If on skin, skin irritation, get

medical advice/attention.

Ingestion: Rinse mouth. Get medical advice/attention.

Inhalation: Remove victim to fresh air and keep at rest in a position comfortable for breathing.

If symptom concerning breath goes out, call a POISON CENTER or doctor.

4.2 Most important symptoms and effects, both acute and delayed

No further relevant information available.

4.3 Indication of any immediate medical attention and special treatment needed

No further relevant information available.

SECTION 5. Fire-fighting measures

5.1 Extinguishing Media:

CO₂, Dry chemical, Foam, Dry sand

5.2 Special hazards arising from the substance or mixture:

No further relevant information available.

5.3 Advice for firefighters:

Wear fire protective cloth and self-contained breathing apparatus, if necessary.

(Contd. on page3)

Revision date: April 1, 2020



page 3/6 BeautiSealant PASTE

Printing date: June 10, 2020

(Contd. of page2)

SECTION 6. Accidental release measures

6.1 Personal precautions, protective equipment and emergency procedures:

Avoid eye and skin contact.

6.2 Environmental Precautions:

Send to approved treatment/disposal company or dispose according to local, state and federal regulations.

6.3 Methods and material for containment and cleaning up:

Wipe up and discard in a suitable container.

6.4 Reference to other section:

See Section 7 for information on safe handling.

See Section 8 for information on personal protection equipment.

See Section 13 for disposal information.

SECTION 7. Handling and storage

7.1 Precautions for safe handling:

Avoid inhalation of grindings and prolonged skin contact with uncured resin.

7.2 Conditions for safe storage, including any incompatibilities:

Store in a cool and dark area.

7.3 Specific end use(s):

No further relevant information available.

SECTION 8. Exposure controls/personal protection

8.1 Control parameters:

Ingredients with limit values that require monitoring at the workplace:

The product does not contain any relevant quantities of materials with critical values that have to be monitored at the workplace.

8.2 Exposure controls:

Respiratory Protection:

Not required (use dust mask while grinding or finishing.)

Skin Protection: Hand Protection

The glove material has to be impermeable and resistant to the product/ the

substance/ the preparation.

Due to missing tests no recommendation to the glove material can be given

for the product/ the preparation/the chemical mixture.

Selection of the glove material on consideration of the penetration times,

rates of diffusion and the degradation.

Material of gloves

The selection of the suitable gloves does not only depend on the material, but also on further marks of quality and varies from manufacturer to manufacturer. As the product is a preparation of several substances, the resistance of the glove material can not be calculated in advance and has therefore to be checked

prior to the application.

(Contd. on page4)

Version Number 5 Revision date: April 1, 2020



page 4/6 **BeautiSealant PASTE**

(Contd. of page3)

Printing date: June 10, 2020

· Penetration time of glove material

The exact break through time has to be found out by the manufacturer of the protective gloves and has to be observed.

· For the permanent contact of a maximum of 15 minutes gloves made of the following materials are suitable:

> Butyl rubber, BR Nitrile rubber, NBR

Eye Protection: Safety goggles

SECTION 9. Physical and chemical properties

Information on basic physical and chemical properties

Appearance/Odor/Colour: Low viscosity paste of the lemon yellow.

Odour threshold Not determined. Not determined. pΗ Melting point/freezing point Not determined. **Boiling Point:** Not determined. Flash point: Not determined. Evaporation rate Not determined. Flammability (solid, gas) Not applicable. Upper/lower flammability or explosive limits Not determined. Not determined. Vapour pressure Vapour density Not determined. Relative Density: ≈ 1.2 (water=1)

Solubility: water solubility Insoluble

Partition coefficient: n-octanol/water Not determined. Auto-ignition temperature Not determined. Decomposition temperature Not determined. Viscosity Not determined. Explosive properties Not applicable. Oxidising properties Not applicable.

9.2 Other information

No further relevant information available.

SECTION 10. Stability and reactivity

10.1 Reactivity:

Polymerization will occur with heat, light.

10.2 Chemical stability:

Stable under normal temperatures and pressures.

10.3 Possibility of hazardous reactions:

No dangerous reactions known.

10.4 Condition to Avoid:

Direct sunlight, high temperature and heating.

10.5 Incompatible materials:

Strong oxidizing materials.

(Contd. on page5)

Version Number 5 Revision date: April 1, 2020



page 5/6
BeautiSealant
PASTE

Printing date: June 10, 2020

(Contd. of page4)

10.6 Hazardous Decomposition Products:

None under normal conditions of storage and use.

SECTION 11. Toxicological information

11.1 Information on toxicological effects:

Acute toxicity: Based on available data, the classification criteria are not met.

Skin corrosion/irritation: Skin Irrit. 2; H315 Causes skin irritation.

Eye damage/irritation: Eye Irrit. 2; H319 Causes serious eye irritation.

Sensitization to the respiratory tract:

Based on available data, the classification criteria are not met.

Skin sensitization: Skin Sens. 1; H317 May cause an allergic skin reaction.

Germ cell mutagenicity/Genotoxicity:

Based on available data, the classification criteria are not met.

Carcinogenicity: Based on available data, the classification criteria are not met. Reproductive toxicity: Based on available data, the classification criteria are not met.

Effects on or via lactation: Lack of data.

Specific target organ toxicity (single exposure):

Based on available data, the classification criteria are not met.

Specific target organ toxicity (repeated exposure):

Based on available data, the classification criteria are not met.

Aspiration hazard: Based on available data, the classification criteria are not met.

SECTION 12. Ecological information

12.1 Toxicity:

No further relevant information available.

12.2 Persistence and degradability:

No further relevant information available.

12.3 Bioaccumulative potential:

No further relevant information available.

12.4 Mobility in soil:

No further relevant information available.

12.5 Results of PBT and vPvB assessment:

Not applicable.

12.6 Other adverse effects:

No further relevant information available.

SECTION 13. Disposal considerations

13.1 Waste treatment methods:

Dispose of contents/container to in accordance with local/regional/national/international regulations.

SECTION 14. Transport information

14.1 UN number: Void14.2 UN proper shipping name: Void

(Contd. on page6)

Revision date: April 1, 2020



page 6/6 BeautiSealant PASTE

Printing date: June 10, 2020

(Contd. of page5)

14.3 Transport hazard class(es): Void14.4 Packing group: Void

14.5 Environmental hazards: No further relevant information available.

14.6 Special precautions for user: Not applicable.

14.7 Transport in bulk according to Annex II of MARPOL73/78 and the IBC Code:

Not applicable.

SECTION 15. Regulatory information

15.1 Safety, health and environmental regulations/legislation specific for the substance or mixture:

• EU REGULATIONS: See Section 2

· Other regulations, limitations and prohibitive regulations:

The product is a medical device according to the EC-directive 93/42/EEC.

15.2 Chemical safety assessment:

A Chemical Safety Assessment has not been carried out.

SECTION 16. Other information

This product is intended for use by dental professionals. (instrument/material)

Relevant phrases:

H315 Causes skin irritation.

H317 May cause an allergic skin reaction.

H319 Causes serious eye irritation.

Abbreviations and acronyms:

EINESC: European Inventory of Existing Commercial Chemical Substances CAS: Chemical Abstracts Service (division of the American Chemical Society)

PBT: Persistent, Bioaccumulative and Toxic vPvB: very Persistent and very Bioaccumulative

Version Number 5 Revision date: April 1, 2020