

Safety data sheet

according to 1907/2006/EC, Article 31

Printing date 23.07.2019

Version number 1

Revision: 23.07.2019

SECTION 1: Identification of the substance/mixture and of the company/ undertaking

- 1.1 Product identifier
- Trade name: **POUDRE AIR N GO CLASSIC NEUTRE / AIR N GO CLASSIC NEUTRAL POWDER**
- 1.2 Relevant identified uses of the substance or mixture and uses advised against
No further relevant information available.
- Application of the substance / the preparation : Powder for dental polishing
- 1.3 Details of the supplier of the safety data sheet
- Manufacturer/Supplier:
SATELEC
Z.I. du Phare – 17 avenue Gustave Eiffel
33700 MERIGNAC - FRANCE
Tel: + 33 (0)5 56.34.06.07
Fax :+ 33 (0)5 56.34.92.92
Email :olivier.petit@acteongroup.com
- 1.4 Emergency telephone number:
France (INRS 24h/24) - Tel : +33 (0)1 45 42 59 59
Ireland - Tel : 00 353 1 8092568 - 00 353 1 8379964 (24h/24) EU Tel: 112
Belgique - Tel : 32 070/245 245
UK : Call NHS 999 in a medical emergency. If you urgently need medical help or advice but it's not a life-threatening situation, call 111 (24/7). The NHS 111 number is currently only available in certain areas. If you're outside of these areas, you should call NHS Direct on 0845 4647

SECTION 2: Hazards identification

- 2.1 Classification of the substance or mixture
- Classification according to Regulation (EC) No 1272/2008 :
The product is not classified, according to the CLP regulation.
- 2.2 Label elements
- Labelling according to Regulation (EC) No 1272/2008
The product is classified and labelled according to the CLP regulation.
- Hazard pictograms Not classified as hazardous
- Signal word Not classified as hazardous
- Hazard statements Not classified as hazardous
- Precautionary statements
P102 Keep out of reach of children.
P280 Wear protective gloves / eye protection / face protection.
P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
P337+P313 If eye irritation persists: Get medical advice/attention.
- 2.3 Other hazards
- Results of PBT and vPvB assessment
- PBT: Not applicable.
- vPvB: Not applicable.

SECTION 3: Composition/information on ingredients

- 3.2 Chemical characterisation: Mixtures
- Description: Mixture: consisting of the following components.

(Contd. on page 2)

— GB —

Safety data sheet

according to 1907/2006/EC, Article 31

Printing date 23.07.2019

Version number 1

Revision: 23.07.2019

Trade name: **POUDRE AIR N GO CLASSIC NEUTRE / AIR N GO CLASSIC NEUTRAL POWDER**

(Contd. of page 1)

| | | |
|---|--------------------------|---------|
| · Dangerous components: | | |
| CAS: 144-55-8 EINECS: 205-633-8 Reg.nr.: 01-2119457606-32 | sodium hydrogencarbonate | 50-100% |

· SVHC No

· Additional information: For the wording of the listed risk phrases, refer to section 16.

SECTION 4: First aid measures

- 4.1 Description of first aid measures
- General information: No special measures required.
- After inhalation:
Supply fresh air. If required, provide artificial respiration. Keep patient warm. Consult doctor if symptoms persist.
In case of unconsciousness place patient in a stable laying down side position for transportation.
- After skin contact:
Immediately wash with water and soap and rinse thoroughly.
Generally the product does not irritate the skin.
- After eye contact: Rinse opened eye for several minutes under running water.
- After swallowing:
Do not induce vomiting; call for medical help immediately.
If symptoms persist consult a doctor.
- 4.2 Most important symptoms and effects, both acute and delayed :
No further relevant information available.
- 4.3 Indication of any immediate medical attention and special treatment needed :
No further relevant information available.

SECTION 5: Firefighting measures

- 5.1 Suitable extinguishing agents:
CO₂, powder or water spray. Fight larger fires with water spray or alcohol resistant foam.
- 5.2 Special hazards arising from the substance or mixture : No further relevant information available.
- 5.3 Advice for firefighters
- Protective equipment: No special measures required.

SECTION 6: Accidental release measures

- 6.1 Personal precautions, protective equipment and emergency procedures :
Wear a respiratory protective device.
Remove persons from the danger area.
- 6.2 Environmental precautions: Do not allow the product to enter sewers/ surface water or ground water.
- 6.3 Methods and material for containment and cleaning up: Pick up mechanically.
- 6.4 Reference to other sections :
See Section 7 for information on safe handling.
See Section 8 for information on personal protective equipment.
See Section 13 for disposal information.

SECTION 7: Handling and storage

- 7.1 Precautions for safe handling :
Prevent formation of dust.

(Contd. on page 3)

GB

Safety data sheet according to 1907/2006/EC, Article 31

Printing date 23.07.2019

Version number 1

Revision: 23.07.2019

Trade name: **POUDRE AIR N GO CLASSIC NEUTRE / AIR N GO CLASSIC NEUTRAL
POWDER**

(Contd. of page 2)

- Keep receptacles tightly sealed.
- Ensure good ventilation/exhaustion at the workplace.
- Information about fire - and explosion protection: No special measures required.
- 7.2 Conditions for safe storage, including any incompatibilities
- Storage:
- Requirements to be met by storerooms and receptacles: Store only in the original receptacle.
- Information about storage in one common storage facility: Not required.
- Further information about storage conditions: Protect from heat and direct sunlight.
- Recommended storage temperature: 5°C-25°C
- 7.3 Specific end use(s) : No further relevant information available.

SECTION 8: Exposure controls/personal protection

- Additional information about design of technical facilities: No further data; see section 7.
- 8.1 Control parameters
- Ingredients with limit values that require monitoring at the workplace:
The product does not contain any significant quantities of materials with critical values that have to be monitored at the workplace.
- Additional information:
The most current valid lists have been used as a basis for the production of this document.
- 8.2 Exposure controls
- Personal protective equipment:
- General protective and hygienic measures:
Do not eat and drink while not using
The usual precautionary measures are to be adhered to when handling chemicals.
- Respiratory protection:



In case of brief exposure or low pollution use respiratory filter device. In case of intensive or longer exposure use self-contained respiratory protective device.

- Protection of hands:



Protective gloves

- Gloves material : Nitrile rubber, NBR
- Eye protection:



Safety glasses

SECTION 9: Physical and chemical properties

- 9.1 Information on basic physical and chemical properties
- General Information
- Appearance:
- Form: Powder
- Colour: White

(Contd. on page 4)

Safety data sheet

according to 1907/2006/EC, Article 31

Printing date 23.07.2019

Version number 1

Revision: 23.07.2019

Trade name: **POUDRE AIR N GO CLASSIC NEUTRE / AIR N GO CLASSIC NEUTRAL POWDER**

(Contd. of page 3)

| | |
|--|---|
| · Odour: | Odourless |
| · Odour threshold: | Not determined. |
| · pH-value: | Not applicable. |
| · Change in condition | |
| Melting point/freezing point: | Undetermined. |
| Initial boiling point and boiling range: | Undetermined. |
| · Flash point: | Not applicable. |
| · Flammability (solid, gas): | Not determined. |
| · Ignition temperature: | Not applicable. |
| · Decomposition temperature: | Not determined. |
| · Auto-ignition temperature: | Product is not self igniting. |
| · Explosive properties: | Product does not present an explosion hazard. |
| · Explosion limits: | |
| Lower: | Not determined. |
| Upper: | Not determined. |
| · Vapour pressure: | Not applicable. |
| · Density: | Not determined. |
| · Relative density : | Not determined. |
| · Vapour density : | Not applicable. |
| · Evaporation rate : | Not applicable. |
| · Solubility in / Miscibility with water: | Fully soluble. |
| · Segregation coefficient (n-octanol/water): | Not determined. |
| · Viscosity: | |
| Dynamic: | Not applicable. |
| Kinematic: | Not applicable. |
| Kinematic: | Not applicable. |
| · Solvent content: | |
| COV (CH) : | 0.00 % |
| · 9.2 Other information : | No further relevant information available. |

SECTION 10: Stability and reactivity

- 10.1 Reactivity No further relevant information available.
- 10.2 Chemical stability
- Thermal decomposition / conditions to be avoided: No decomposition if used according to specifications.
- 10.3 Possibility of hazardous reactions : No dangerous reactions known.
- 10.4 Conditions to avoid : No further relevant information available.
- 10.5 Incompatible materials: No further relevant information available.
- 10.6 Hazardous decomposition products: No dangerous decomposition products known.

(Contd. on page 5)

Safety data sheet
according to 1907/2006/EC, Article 31

Printing date 23.07.2019

Version number 1

Revision: 23.07.2019

Trade name: **POUDRE AIR N GO CLASSIC NEUTRE / AIR N GO CLASSIC NEUTRAL POWDER**

(Contd. of page 4)

SECTION 11: Toxicological information

- 11.1 Information on toxicological effects
- Acute toxicity Based on available data, the classification criteria are not met.

· LD/LC50 values relevant for classification:

CAS: 144-55-8 sodium hydrogencarbonate

| | | |
|------|------|-------------------|
| Oral | LD50 | 4,220 mg/kg (rat) |
|------|------|-------------------|

- Primary irritant effect:
- On the skin: Based on available data, the classification criteria are not met.
- On the eyes: Based on available data, the classification criteria are not met.
- Respiratory or skin sensitization Based on available data, the classification criteria are not met.
- CMR effects (carcinogenicity, mutagenicity and toxicity for reproduction)
- Germ cell mutagenicity Based on available data, the classification criteria are not met.
- Carcinogenicity Based on available data, the classification criteria are not met.
- Reproductive toxicity Based on available data, the classification criteria are not met.
- STOT-single exposure Based on available data, the classification criteria are not met.
- STOT-repeated exposure Based on available data, the classification criteria are not met.
- Aspiration hazard Based on available data, the classification criteria are not met.

SECTION 12: Ecological information

- 12.1 Toxicity
- Aquatic toxicity: No further relevant information available.
- 12.2 Persistence and degradability : No further relevant information available.
- 12.3 Bioaccumulative potential : No further relevant information available.
- 12.4 Mobility in soil : No further relevant information available.
- Additional ecological information:
- General notes:
Water hazard class 1 (German Regulation) (Self-assessment): slightly hazardous for water.
Do not allow undiluted product or large quantities of it to reach ground water, water course or sewage system.
- 12.5 Results of PBT and vPvB assessment
- PBT: Not applicable.
- vPvB: Not applicable.
- 12.6 Other adverse effects : No further relevant information available.

SECTION 13: Disposal considerations

- 13.1 Waste treatment methods
- Recommendation :



Must not be disposed of together with household garbage. Do not allow product to reach sewage system.

- Uncleaned packaging:
- Recommendation: Disposal must be made according to official regulations.
- Recommended cleansing agents: Water, if necessary together with cleansing agents.

(Contd. on page 6)

Safety data sheet

according to 1907/2006/EC, Article 31

Printing date 23.07.2019

Version number 1

Revision: 23.07.2019

Trade name: **POUDRE AIR N GO CLASSIC NEUTRE / AIR N GO CLASSIC NEUTRAL POWDER**

(Contd. of page 5)

SECTION 14: Transport information

| | |
|---|---|
| · 14.1 UN-Number · ADR, ADN, IMDG, IATA | Not classified as hazardous for transport |
| · 14.2 UN proper shipping name · ADR, ADN, IMDG, IATA | Not classified as hazardous for transport |
| · 14.3 Transport hazard class(es) · ADR, ADN, IMDG, IATA · Class | Not classified as hazardous for transport |
| · 14.4 Packing group · ADR, IMDG, IATA | Not classified as hazardous for transport |
| · 14.5 Environmental hazards: | Not applicable. |
| · 14.6 Special precautions for user | Not applicable. |
| · 14.7 Transport in bulk according to Annex II of Marpol and the IBC Code | Not applicable. |
| · UN "Model Regulation": | Not classified as hazardous for transport |

SECTION 15: Regulatory information

- 15.1 Safety, health and environmental regulations/legislation specific for the substance or mixture
- Directive 2012/18/EU
- Named dangerous substances - ANNEX I None of the ingredients is listed.
- 15.2 Chemical safety assessment: A Chemical Safety Assessment has not been carried out.

SECTION 16: Other information

This information is based on our present knowledge. However, this shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.

| |
|---|
| · Classification according to Regulation (EC) No 1272/2008 |
| Self-reactive substances and mixtures On basis of test data |

· Abbreviations and acronyms:

ADR: Accord européen sur le transport des marchandises dangereuses par Route (European Agreement concerning the International Carriage of Dangerous Goods by Road)
 IMDG: International Maritime Code for Dangerous Goods
 IATA: International Air Transport Association
 GHS: Globally Harmonised System of Classification and Labelling of Chemicals
 EINECS: European Inventory of Existing Commercial Chemical Substances
 ELINCS: European List of Notified Chemical Substances
 CAS: Chemical Abstracts Service (division of the American Chemical Society)
 VOCV: Lenkungsabgabe auf flüchtigen organischen Verbindungen, Schweiz (Swiss Ordinance on volatile organic compounds)
 LC50: Lethal concentration, 50 percent
 LD50: Lethal dose, 50 percent
 PBT: Persistent, Bioaccumulative and Toxic
 SVHC: Substances of Very High Concern
 vPvB: very Persistent and very Bioaccumulative