

# Safety data sheet

according to Regulation (EC) No 1907/2006

Trade name: **b.clean**

Issue/Revision: 02.04.2025

Version: 1.0

Replaces version: -

## SECTION 1: Identification of the substance/mixture and of the company/undertaking

### 1.1. Product identifier

Trade name: **b.clean**  
UFI: not applicable

### 1.2. Relevant identified uses of the substance or mixture and uses advised against

Relevant identified uses: Cleaning and disinfecting agent  
Intended purpose: Solution ready for use for the decontamination and prevention of biofilm formation in procedural water lines of medical and dental units.  
Uses advised against: None at intended use.  
Note: The product is intended for professional users.

### 1.3. Details of the supplier of the safety data sheet

Supplier: NSK Europe GmbH  
Elly-Beinhorn-Straße 8  
D-65760 Eschborn (Germany)  
Telephone: +49 6196 77606-0  
Telefax: +49 6196 77606-29  
E-Mail: info@nsk-europe.de  
  
Manufacturer: ALPRO MEDICAL GMBH  
Mooswiesenstraße 9  
D-78112 St. Georgen (Germany)  
Telephone: +49 7725 9392-0  
Telefax: +49 7725 9392-91  
E-mail: info@alpro-medical.de  
Internet: www.alpro-medical.de

E-mail address for the competent person responsible for the safety data sheet: doku@alpro-medical.de

### 1.4. Emergency telephone number

In-house emergency telephone number: +49 7725 9392-0  
Monday – Friday from 08:00 am to 04:30 pm (UTC+1);  
for chemical information and legal information on hazardous substances only  
  
Poison centre: +49 761 19240  
Poisoning information centre, Freiburg, Germany  
(24 h / 7 d), English is spoken  
  
National Poisons Information Service (UK): +44 344 892 0111  
National Poisons Information Service (NPIS)  
(24 h / 7 d), Medical Professionals Only

## SECTION 2: Hazards identification

### 2.1. Classification of the substance or mixture

Classification in accordance with Regulation (EC) No 1272/2008 [CLP]

The mixture is not classified as hazardous in the sense of Regulation (EC) No 1272/2008 [CLP].

# Safety data sheet

according to Regulation (EC) No 1907/2006

Trade name: **b.clean**

Issue/Revision: 02.04.2025

Version: 1.0

Replaces version: -

## 2.2. Label elements

### Label elements in accordance with Regulation (EC) No 1272/2008 [CLP]

The mixture does not require labelling in the sense of Regulation (EC) No 1272/2008 [CLP].

Hazard pictograms: -

Signal word: -

Hazard components  
for labelling: -

H-phrases: -

P-phrases: -

## 2.3. Other hazards

The substances in the mixture do not meet the PBT/vPvB criteria according to Regulation (EC) No 1907/2006, Annex XIII (see SECTION 12.5.)

The substances in the mixture have no endocrine disrupting properties according to Regulation (EC) No. 1907/2006, Annex XIV (see SECTION 11 and SECTION 12.6.). They are not on the list of substances of very high concern for authorization according to Regulation (EC) No. 1907/2006, Article 59, Paragraph 10).

No further hazards known.

## SECTION 3: Composition/information on ingredients

### 3.1. Substances

Not applicable.

### 3.2. Mixtures

Chemical characterisation: Mixture of propylene glycol, parabens, biguanides and complexing agents in aqueous solution.

#### Hazardous ingredients

Does not contain hazardous ingredients above concentration limit values fixed in Regulation (EC) No 1907/2006.

## SECTION 4: First aid measures

### 4.1. Description of first aid measures

General information: First aider: Pay attention to self-protection!

Following inhalation: Supply fresh air. In case of complaints consult a physician.

Following skin contact: No special measures necessary.

Following eye contact: Flush eyes with flowing water for several minutes holding eyelids apart.  
Remove contact lenses, if present and easy to do.

Following ingestion: Rinse mouth with water. Let drink water.

### 4.2. Most important symptoms and effects, both acute and delayed

None known.

### 4.3. Indication of any immediate medical attention and special treatment needed

The product contains a source of phenylalanine.

# Safety data sheet

according to Regulation (EC) No 1907/2006

Trade name: **b.clean**

Issue/Revision: 02.04.2025

Version: 1.0

Replaces version: -

## SECTION 5: Firefighting measures

### 5.1. Extinguishing media

Suitable extinguishing media: Water spray jet, alcohol resistant foam, extinguishing powder, carbon dioxide (CO<sub>2</sub>)

Unsuitable extinguishing media: Full water jet

### 5.2. Special hazards arising from the substance or mixture

Hazardous combustion products: Carbon monoxide (CO), carbon dioxide (CO<sub>2</sub>)

### 5.3. Advice for firefighters

Special protective equipment: Not necessary

Further information: Cool endangered containers with water spray jet.

## SECTION 6: Accidental release measures

### 6.1. Personal precautions, protective equipment and emergency procedures

#### For non-emergency personnel

Use personal protective equipment. See SECTION 8.2.

Special danger of slipping by leaked/spilled product.

#### For emergency responders

Use personal protective equipment. See SECTION 8.2.

### 6.2. Environmental precautions

Not necessary

### 6.3. Methods and material for containment and cleaning up

#### Containment

Not necessary

#### Cleaning up

Wipe up small amounts with absorbent material (e.g. cloth, fleece). Absorb large amounts with liquid-binding material (sand, diatomaceous earth, universal binder, sawdust).

#### Other information

Inappropriate containment and cleaning methods are not known.

### 6.4. Reference to other sections

Information on safe handling see SECTION 7.1.

Information on personal protective equipment see SECTION 8.2.

Information on disposal see SECTION 13.

## SECTION 7: Handling and storage

### 7.1. Precautions for safe handling

#### Precautions

Not necessary

#### Advice on general occupational hygiene

When using do not eat, drink or smoke. Wash hands before breaks and at end of work. Keep away from medicines, food, feed, cosmetics and stimulants.

# Safety data sheet

according to Regulation (EC) No 1907/2006

Trade name: **b.clean**

Issue/Revision: 02.04.2025

Version: 1.0

Replaces version: -

## 7.2. Conditions for safe storage, including any incompatibilities

Requirements for storage rooms and vessels: Keep only in the original container.

Advice on common storage: Not necessary

Further information on storage conditions: Do not store above 30°C.

## 7.3. Specific end use(s)

Apart from the uses mentioned in section 1.2 no other specific end uses are stipulated.

### Industry and sector specific guidance

[DE] TRGS 525 – Hazardous substances in medical care facilities (Section 7 Activities with disinfectants); Issue: September 2014; Source: GMBI 2014 page 1294-1307 of 13.10.2014 [No 63]; 10.07.2015 [No 27]; [www.baua.de](http://www.baua.de)

[DE] DGUV rules 207-206 – Prevention of chemical risks when handling disinfectants in health service, Issue: 2016.12; Source: [www.https://dguv.de/publikationen](https://dguv.de/publikationen)

## SECTION 8: Exposure controls/personal protection

### 8.1. Control parameters

#### Occupational exposure limits

Does not contain substances above concentration limits fixing an occupational exposure limit.

#### Biological limit values

Does not contain substances above concentration limits fixing a biological limit value.

#### Information on monitoring procedures

BS EN 482:2021-04; Title: Workplace exposure - Procedures for the determination of the concentration of chemical agents;  
British version of EN 482:2021

BS EN 689:2018; Title: Workplace atmospheres - Guidance for the assessment of exposure by inhalation to chemical agents for comparison with limit values and measurement strategy;  
British version of EN 689:2018

BS EN 14042:2003-04-24; Title: Workplace atmospheres - Guide for the application and use of procedures for the assessment of exposure to chemical and biological agents;  
British version of EN 14042:2003)

CEN/TR 17055:2017; Title: Workplace exposure. Measurement of chemical agents complying with the requirements given in EN 482 and either one of EN 838, EN 1076, EN 13205, EN 13890 and EN 13936.  
Choice of procedures

ISO TR 14294:2011; Title: Workplace atmospheres. Measurement of dermal Exposure. Principles and methods

### 8.2. Exposure controls

#### Appropriate engineering controls

##### Technical and organisational protective measures

No special protective measures necessary.

#### Personal protective equipment

Eye/face protection: Not necessary when used as intended.

Skin protection:

Hand protection: Not necessary when used as intended.

# Safety data sheet

according to Regulation (EC) No 1907/2006

Trade name: **b.clean**

Issue/Revision: 02.04.2025

Version: 1.0

Replaces version: -

Other skin protection: Not necessary when used as intended.

Respiratory protection: Not necessary when used as intended.

Thermal hazards: No special protective measures necessary.

## Environmental exposure controls

Not necessary

## SECTION 9: Physical and chemical properties

Appearance/physical state:	clear, turquoise liquid	
Odour:	apple	
Odour threshold:	no data available	
Melting point/freezing point:	no data available	
Initial boiling point and boiling range:	no data available	
Flammability:	not applicable	
Lower explosive limit:	not applicable	
Upper explosive limit:	not applicable	
Flash point:	not applicable	
Auto-ignition temperature:	no data available	
pH (undiluted):	5.5 – 6.5	(20 °C)
Kinematic viscosity:	no data available	
Solubility in water:	completely soluble	
Partition coefficient: n-octanol/water	not applicable	
Vapour pressure:	no data available	(... °C)
Density:	0.995 – 1.005 g/cm <sup>3</sup>	(20 °C)
Relative vapour density:	no data available	
Particle characteristics:	not applicable	

### 9.2. Other information

#### 9.2.1. Information with regard to physical hazard classes

Explosive substances/mixtures and products containing explosives:	not applicable
Flammable gases:	not applicable
Aerosols:	not applicable
Oxidising gases:	not applicable
Gases under pressure:	not applicable
Flammable liquids:	not applicable
Flammable solids:	not applicable
Self-reactive substances and mixtures:	not applicable
Pyrophoric liquids:	not applicable
Pyrophoric solids:	not applicable

# Safety data sheet

according to Regulation (EC) No 1907/2006

Trade name: **b.clean**

Issue/Revision: 02.04.2025

Version: 1.0

Replaces version: -

Self-heating substances and mixtures:	not applicable
Substances and mixtures, which emit flammable gases in contact with water:	not applicable
Oxidising liquids:	not applicable
Oxidising solids:	not applicable
Organic peroxides:	not applicable
Substances and mixtures corrosive to metals:	not applicable
Desensitised explosives:	not applicable

## 9.2.2. Other safety characteristics

Electrical conductivity (undiluted): 920-1020  $\mu\text{S}/\text{cm}$  (20 °C)

## SECTION 10: Stability and reactivity

### 10.1. Reactivity

No hazardous reactions when handled and stored as intended.

### 10.2. Chemical stability

The product is stable when handled and stored as intended.

### 10.3. Possibility of hazardous reactions

None known

### 10.4. Conditions to avoid

None known

### 10.5. Incompatible materials

No material incompatibilities when used as intended.

### 10.6. Hazardous decomposition products

Does not decompose when used as intended.

## SECTION 11: Toxicological information

### 11.1. Information on toxicological effects

#### Acute toxicity

Product

Acute toxicity - oral:	Acute Toxicity Estimate $\text{ATE}_{\text{mix}} > 2000 \text{ mg/kg}$ => no classification
Acute toxicity - dermal:	Acute Toxicity Estimate $\text{ATE}_{\text{mix}} > 2000 \text{ mg/kg}$ => no classification
Acute toxicity - inhalation:	Acute Toxicity Estimate $\text{ATE}_{\text{mix}} > 20 \text{ mg/l}$ => no classification

#### Skin corrosion/irritation

Product

No classification. [calculation method]

# Safety data sheet

## according to Regulation (EC) No 1907/2006

Trade name: **b.clean**  
Issue/Revision: 02.04.2025

Version: 1.0  
Replaces version: -

### Serious eye damage/irritation

Product

No classification. [calculation method]

### Respiratory or skin sensitisation

Product

No classification. [calculation method]

### Germ cell mutagenicity

Product

No data available.

### Carcinogenicity

Product

No data available.

### Reproductive toxicity

Product

No data available.

### STOT-single exposure

Product

No classification. [calculation method]

### STOT-repeated exposure

Product

No data available.

### Aspiration hazard

Product

No data available.

### Information on other hazards

Endocrine disrupting properties:

No substances are contained that have endocrine disrupting properties for humans.

## SECTION 12: Ecological information

### 12.1. Toxicity

No classification. [calculation method]

### 12.2. Persistence and degradability

Biodegradability:

The product is biodegradable according to OECD criteria. The statement has been derived from the properties of the ingredients.

### 12.3. Bioaccumulative potential

No data available.

### 12.4. Mobility in soil

No data available.

# Safety data sheet

according to Regulation (EC) No 1907/2006

Trade name: **b.clean**

Issue/Revision: 02.04.2025

Version: 1.0

Replaces version: -

## 12.5. Results of PBT and vPvB assessment

The substances in the mixture do not meet the PBT/vPvB criteria according to Regulation (EC) No 1907/2006, Annex XIII.

## 12.6. Endocrine disrupting properties

No substances are contained that have endocrine disrupting properties for non-target organisms.

## 12.7. Other adverse effects

No data available.

## SECTION 13: Disposal considerations

### 13.1. Waste treatment methods

#### Disposal of the product

Product residues must be disposed of as non-hazardous waste in compliance with the Directive 2008/98/EC on waste as well as national and regional regulations. Do not mix with other waste materials.

Waste codes / waste designations according to EWC

Product residues: 16 10 02 aqueous liquid wastes other than those mentioned in 16 10 01

#### Disposal of the packaging

Waste codes / waste designations according to EWC

Contaminated packaging: 15 01 02 plastic packaging

Recommendation

Contaminated packaging must be emptied optimally and can be recycled.

## SECTION 14: Transport information

### 14.0. Transport classification

No dangerous good in sense of the transport regulations in road traffic (ADR), railway traffic (RID), inland waterway traffic (ADN), maritime traffic (IMDG-Code) and air traffic (ICAO-TI/IATA-DGR).

### 14.1. UN number

-

### 14.2. UN proper shipping name

-

### 14.3. Transport hazard class(es)

-

### 14.4. Packing group

-

### 14.5. Environmental hazards

-

### 14.6. Special precautions for user

-



# Safety data sheet

according to Regulation (EC) No 1907/2006

Trade name: **b.clean**

Issue/Revision: 02.04.2025

Version: 1.0

Replaces version: -

## 14.7. Maritime transport in bulk according to IMO instruments

-

## 14.8. Further information

-

## SECTION 15: Regulatory information

### 15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture

#### EU-Regulations

REGULATION (EC) No 1005/2009 on substances that deplete the ozone layer

not applicable

REGULATION (EC) No 2019/1021 on persistent organic pollutants

not applicable

REGULATION (EU) No 649/2012 concerning the export and import of hazardous chemicals

not applicable

REGULATION (EU) No 648/2004 on detergents

Disinfectants

Perfumes

DIRECTIVE 2012/18/EU (Seveso III Directive) on the control of major-accident hazards involving dangerous substances, amending and subsequently repealing Council Directive 96/82/EC

not applicable

DIRECTIVE 2010/75/EU on industrial emissions (integrated pollution prevention and control)

VOC content: < 2 %

REACH – List of substances subject to authorisation (Annex XIV)

not applicable

REACH – Restrictions on the manufacture, placing on the market and use of certain dangerous substances, mixtures and articles (Annex XVII)

not applicable

COUNCIL DIRECTIVE 94/33/EC on the protection of young people at work

not applicable

COUNCIL DIRECTIVE 92/85/EEC on the introduction of measures to encourage improvements in the safety and health at work of pregnant workers and workers who have recently given birth or are breastfeeding

not applicable

### 15.2. Chemical safety assessment

For this mixture no chemical safety assessment has been carried out.

## SECTION 16: Other information

### 16.1. Full text of hazard classes and H-phrases

#### Hazard classes

-

# Safety data sheet

according to Regulation (EC) No 1907/2006

Trade name: **b.clean**

Issue/Revision: 02.04.2025

Version: 1.0

Replaces version: -

## H-phrases (Hazard statements)

-

## 16.2. Abbreviations and acronyms

ADN	<u>A</u> ccord européen relatif au transport international des marchandises <u>d</u> angereuses par voie de <u>n</u> avigation intérieure (European Agreement concerning the International Carriage of Dangerous Goods by Inland Waterways)
ADR	<u>A</u> ccord européen relatif au transport international des marchandises <u>d</u> angereuses par <u>r</u> oute (European Agreement concerning the International Carriage of Dangerous Goods by Road)
BS	<u>B</u> ritish <u>S</u> tandards
CAS	<u>C</u> hemical <u>A</u> bstracts <u>S</u> ervice
CLP	Regulation on <u>C</u> lassification, <u>L</u> abelling and <u>P</u> ackaging of Substances and Mixtures
[DE]	National German regulations
DGUV	<u>D</u> eutsche <u>G</u> esetzliche <u>U</u> nfall <u>v</u> ersicherung (English: German statutory accident insurance)
EC	<u>E</u> uropean <u>C</u> ommunity
EEC	<u>E</u> uropean <u>E</u> conomic <u>C</u> ommunity
EN	European Standard
EU	<u>E</u> uropean <u>U</u> nion
EWC	<u>E</u> uropean <u>W</u> aste <u>C</u> atalogue
GMBI	<u>G</u> emeinsames <u>M</u> inisterial <u>b</u> latt (English: Joint Ministerial Gazette)
IATA-DGR	<u>I</u> nternational <u>A</u> ir <u>T</u> ransport <u>A</u> ssociation - <u>D</u> angerous <u>G</u> oods <u>R</u> egulations
ICAO-TI	Technical Instructions For The Safe Transport of Dangerous Goods by Air
IMDG-Code	<u>I</u> nternational <u>M</u> aritime <u>C</u> ode for <u>D</u> angerous <u>G</u> oods
OECD	<u>O</u> rganization for <u>E</u> conomic <u>C</u> o-operation and <u>D</u> evelopment
PBT	<u>P</u> ersistent, <u>b</u> ioaccumulative and <u>t</u> oxic
ppm	<u>P</u> arts <u>p</u> er <u>m</u> illion
REACH	<u>R</u> egistration, <u>E</u> valuation, <u>A</u> uthorisation and Restriction of <u>C</u> hemicals
RID	<u>R</u> èglement concernant le transport <u>I</u> nternational ferroviaire de marchandises <u>D</u> angereuses (Regulations Concerning the International Carriage of Dangerous Goods by Rail)
TRGS	<u>T</u> echnische <u>R</u> egeln für <u>G</u> efahr <u>s</u> toffe (English: Technical Rules for Hazardous Substances)
UN	<u>U</u> nited <u>N</u> ations
UTC	Coordinated Universal Time (French: Temps Universel Coordonné)
VOC	Volatile Organic Compounds
vPvB	<u>V</u> ery persistent and <u>v</u> ery <u>b</u> ioaccumulative

## 16.3. Key literature references and sources for data

- Regulation (EC) No 1907/2006 (REACH), Annex II
- European Chemicals Agency (ECHA) – Guidance on the compilation of safety data sheets; Version 4.0 (December 2022); <https://echa.europa.eu/documents>
- Regulation (EC) No 1272/2008 (CLP regulation)
- European Chemicals Agency (ECHA) – Guidance on Labelling and Packaging in accordance with Regulation (EC) No 1272/2008; Version 4.2 (03/2021); <https://echa.europa.eu/documents>
- European Chemicals Agency (ECHA), Registered substances; <https://echa.europa.eu/information-on-chemicals/registered-substances>
- European Chemicals Agency (ECHA), C&L Classification and Labelling Inventory; <https://echa.europa.eu/information-on-chemicals/cl-inventory-database>
- Institute for Occupational Safety and Health of the German Social Accident Insurance (IFA): GESTIS database on hazardous substances and GESTIS - International limit values for chemical agents; <https://www.dguv.de/ifa/index.jsp>
- German Environmental Agency (Umweltbundesamt), Section IV 2.4: Office of Documentation and Information on Substances Hazardous to Waters RIGOLETTO (catalogue of Substances Hazardous to Waters); <https://webriigoletto.uba.de/rigoletto>

# Safety data sheet

according to Regulation (EC) No 1907/2006

Trade name: **b.clean**

Issue/Revision: 02.04.2025

Version: 1.0

Replaces version: -

## **16.4. Methods according to Article 9 of Regulation (EC) No 1272/2008 for the evaluation of information for classification purposes**

Calculation method according to the criteria in Annex I 1272/2008.

Flash point according to EN ISO 2719:2002.

pH value measurement.

Material compatibility and corrosiveness in practical tests

## **16.5. Training advice**

Provide adequate information, instructions and training for users.

## **16.6. Indication of changes**

New safety data sheet. No changes so far.

---

The information given in the safety data sheet only applies to the described product in connection with its intended use. This information is based on the latest state of our knowledge at the time of revision. In particular, it describes our product under the aspect of its hazards and safety measures to be taken. It does not constitute any guarantee of product properties and quality features.

---